



Oral 6 Advanced

The **Onsite Oral 6 Advanced** saliva drug screening test incorporates the latest drug groups from **Australian Standard 4760:2019**. One of the first rapid tests on the Australian market designed to the **NEW 2019 AUSTRALIAN STANDARD**.

Saliva testing is the ideal test choice when 'current impairment' screening is required.

Benefits of Oral 6 Advanced :

- EASY-TO USE AND CONVENIENT
- CLEAR, RAPID RESULTS IN MINUTES!
- COST EFFECTIVE TESTING
- NO TOILET FACILITIES REQUIRED

ONSITE ORAL 6 ADVANCED TESTS FOR THE FOLLOWING:

- AMP 50 ng/mL** (Amphetamine)
- MET 50 ng/mL** (Methamphetamine)
- COC 50 ng/mL** (Cocaine)
- OPI 50 ng/mL** (Opiates)
- THC (Parent) 15 ng/mL** (Marijuana)
- OXY 40 ng/mL** (Oxycodone)

Ease-of-collection

Easy to use , convenient, one step oral saliva test. Simply collect the saliva sample with the test swab , insert the test swab into the test chamber and twist seal the lid. The saliva will run up the test strips.

Independent Test Strips

Each drug has its own test strip which is clearly labelled and colour coded allowing for easy interpretation of the results.

Clear Results

After collection remove the privacy strip and read the test result lines. Each drug group is colour coded for ease of recognition. Results are clear and easy to read and appear within 2-3 minutes.



Includes OXYCODONE and low THC level.

NEW AUSTRALIAN STANDARD 4760:2019 TEST



Australian Designed

Order Code (Box of 25): OD-081



Quality Controls

Test control lines use colloidal gold technology with superior sensitivity, accuracy and precision. Cross-reactivity is pretested against most drugs and other substances.

Test Integrity

Integrity of the saliva sample is maintained due to minimal handling and exposure. The saliva sample is observed at all times preventing adulteration by the test subject.

Test control lines confirm the saliva sample has flowed correctly up the test strips and the test is complete.

Comprehensive instructions and procedure cards are included to ensure testing is conducted correctly.

Result Confirmation


Any presumptive-positive results should be confirmed via LC/MS at a certified NATA laboratory using a second saliva sample.

Australian Designed

Cut-off levels are set for the Australian market reducing discrepancies between initial onsite testing and laboratory confirmation test results.

Onsite Oral 6 Advanced Test Procedure

Allow the test device to reach room temperature [15-30°C (59-86°F)] prior to testing. Do not place anything in the mouth including food, drink, gum, or tobacco products for at least 10 minutes prior to collection of oral fluid specimen.




1. Bring the pouch to room temperature before opening it. Remove the test from the sealed pouch and use it as soon as possible.
2. Remove the test device from the sealed pouch and screw the Collector Cap counterclockwise to pull out the whole piece of collection stick with Sponge from the Collection Chamber. (Step 1)
3. Insert the sponge end of the collection stick into the mouth. Close mouth and gently chew the sponge for saliva excretion. Soak sponge into saliva in mouth and swab the inside of the mouth and tongue to collect oral fluid for up to 3 minutes until the sponge becomes completely soft and fully saturated with saliva. No hard spots should be felt on the sponge when saturated. (Step 2)
4. Remove the sponge from the mouth. With gentle pressure, place the collection stick with saturated sponge into Collection Chamber. (Step 3)


Step 1



Step 2



Step 3




Step 4



5. Screw the Collector Cap clockwise to secure the cap and start the timer. (Step 4)

6. Mark patient ID on the test device. Peel off the label to read test results. Wait for the colour line(s) to appear on the test strips. Read results at 10 minutes. Do not read results after 1 hour. (Step 5)
7. Send a second sample extracted from the collector and placed in a sterilised vial to the laboratory for GC/MS confirmation if necessary.

Step 5



- Negative results—** Indicated by the presence of 2 lines. No matter how dark/light
- Presumptive positive—** only the top 'control' line appears
- Invalid—** the control line fails to appear, usually due to too little oral fluid collected.



Nationally Recognised **HLTPAT005 Collect Specimens for Drugs of Abuse Testing** training is available to ensure onsite drug & alcohol screening personnel are correctly trained to use and administer drug & alcohol test equipment under **Australian Standards AS/NZS4308, AS/NZS4760 & AS3547**. This course has been developed to assist our clients in adhering to their legal responsibilities in the workplace and ensure results stand up in a Court of Law. Available via classroom and distance learning. For further details please see our training page.